



SLC40A1 Antibody, HRP conjugated

Product Code	CSB-PA865089LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP59
Immunogen	Recombinant Human Solute carrier family 40 member 1 protein (148-298AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 40 member 1 (Ferroportin-1) (Iron-regulated transporter 1), SLC40A1, FPN1 IREG1 SLC11A3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC40A1