



ASH1L Antibody, Biotin conjugated

Product Code	CSB-PA865114LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR48
Immunogen	Recombinant Human Histone-lysine N-methyltransferase ASH1L protein (1906-2160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone-lysine N-methyltransferase ASH1L (EC 2.1.1.43) (ASH1-like protein) (huASH1) (Absent small and homeotic disks protein 1 homolog) (Lysine N-methyltransferase 2H), ASH1L, KIAA1420 KMT2H
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ASH1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.