




ugl Antibody

Product Code	CSB-PA865671HA01BGAD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9RC92
Immunogen	Recombinant Bacillus sp. Unsaturated glucuronyl hydrolase protein (1-377AA)
Raised In	Rabbit
Species Reactivity	Bacillus sp.
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Unsaturated glucuronyl hydrolase (UGL) (EC 3.2.1.179) (Glycosaminoglycan hydrolase) (Glycuronidase) (Unsaturated uronic acid hydrolase), ugl
Immunogen Species	Bacillus sp.
Research Area	Others
Target Names	ugl

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ugl antibody at 2.8μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 43 kDa

observed band size: 43 kDa