



TRIM4 Antibody, HRP conjugated

Product Code	CSB-PA866336LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C037
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM4 protein (1-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM4 (EC 2.3.2.27) (RING finger protein 87) (RING-type E3 ubiquitin transferase TRIM4) (Tripartite motif-containing protein 4), TRIM4, RNF87
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.