



TFB2M Antibody, FITC conjugated

Product Code	CSB-PA867154LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5Q4
Immunogen	Recombinant Human Dimethyladenosine transferase 2, mitochondrial protein (30-108AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dimethyladenosine transferase 2, mitochondrial (EC 2.1.1.-) (Hepatitis C virus NS5A-transactivated protein 5) (HCV NS5A-transactivated protein 5) (Mitochondrial 12S rRNA dimethylase 2) (Mitochondrial transcription factor B2) (h-mtTFB) (h-mtTFB2) (hTFB2M) (mtTFB2) (S-adenosylmethionine-6-N', N'-adenosyl(rRNA) dimethyltransferase 2), TFB2M, NS5ATP5
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TFB2M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.