



PIPOX Antibody, Biotin conjugated

Product Code	CSB-PA868381LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0Z9
Immunogen	Recombinant Human Peroxisomal sarcosine oxidase protein (259-337AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal sarcosine oxidase (PSO) (EC 1.5.3.1) (EC 1.5.3.7) (L-pipecolate oxidase) (L-pipecolic acid oxidase), PIPOX, LPIPOX PSO
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIPOX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.