



aroA Antibody, HRP conjugated

Product Code	CSB-PA870830HB01ALAE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9R4E4
Immunogen	Recombinant Agrobacterium sp. 3-phosphoshikimate 1-carboxyvinyltransferase protein (1-455AA)
Raised In	Rabbit
Species Reactivity	Agrobacterium sp.
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-phosphoshikimate 1-carboxyvinyltransferase (EC 2.5.1.19) (5-enolpyruvylshikimate-3-phosphate synthase) (EPSP synthase) (EPSPS), aroA
Immunogen Species	Agrobacterium sp. (strain CP4)
Research Area	Others
Target Names	aroA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.