



aroA Antibody, Biotin conjugated

Product Code	CSB-PA870830HD01ALAE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9R4E4
Immunogen	Recombinant Agrobacterium sp. 3-phosphoshikimate 1-carboxyvinyltransferase protein (1-455AA)
Raised In	Rabbit
Species Reactivity	Agrobacterium sp.
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-phosphoshikimate 1-carboxyvinyltransferase (EC 2.5.1.19) (5-enolpyruvylshikimate-3-phosphate synthase) (EPSP synthase) (EPSPS), aroA
Immunogen Species	Agrobacterium sp. (strain CP4)
Research Area	Others
Target Names	aroA