



aroA Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA870830HD01ALAE |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9R4E4 |
| Immunogen | Recombinant Agrobacterium sp. 3-phosphoshikimate 1-carboxyvinyltransferase protein (1-455AA) |
| Raised In | Rabbit |
| Species Reactivity | Agrobacterium sp. |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | 3-phosphoshikimate 1-carboxyvinyltransferase (EC 2.5.1.19) (5-enolpyruvylshikimate-3-phosphate synthase) (EPSP synthase) (EPSPS), aroA |
| Immunogen Species | Agrobacterium sp. (strain CP4) |
| Research Area | Others |
| Target Names | aroA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |