



ROS1 Antibody, HRP conjugated

Product Code	CSB-PA874445LB01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SJQ6
Immunogen	Recombinant Arabidopsis thaliana Protein ROS1 protein (1-320AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ROS1 (EC 3.2.2.-) (DEMETER-like protein 1) (Repressor of silencing 1), ROS1, DML1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	ROS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.