



CHLP Antibody

Product Code	CSB-PA874935LA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9CA67
Immunogen	Recombinant Arabidopsis thaliana Geranylgeranyl diphosphate reductase, chloroplastic protein (44-467AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Geranylgeranyl diphosphate reductase, chloroplastic (EC 1.3.1.83) (Geranylgeranyl reductase), CHLP
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	CHLP