



SEMA6A Antibody, Biotin conjugated

Product Code	CSB-PA875655YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2E6
Immunogen	Recombinant Human Semaphorin-6A protein (755-990AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Semaphorin-6A (Semaphorin VIA) (Sema VIA) (Semaphorin-6A-1) (SEMA6A-1), SEMA6A, KIAA1368 SEMAQ
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SEMA6A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.