



PCNA Antibody, FITC conjugated

Product Code	CSB-PA878702LC01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9M7Q7
Immunogen	Recombinant Arabidopsis thaliana Proliferating cellular nuclear antigen 1 protein (1-263AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proliferating cellular nuclear antigen 1 (PCNA 1), PCNA, PCNA1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	PCNA