



ELOVL4 Antibody, FITC conjugated

Product Code	CSB-PA880925LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZR5
Immunogen	Recombinant Human Elongation of very long chain fatty acids protein 4 protein (270-314AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongation of very long chain fatty acids protein 4 (EC 2.3.1.199) (3-keto acyl-CoA synthase ELOVL4) (ELOVL fatty acid elongase 4) (ELOVL FA elongase 4) (Very long chain 3-ketoacyl-CoA synthase 4) (Very long chain 3-oxoacyl-CoA synthase 4), ELOVL4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ELOVL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.