



NHEJ1 Antibody, Biotin conjugated

Product Code	CSB-PA881005LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9Q4
Immunogen	Recombinant Human Non-homologous end-joining factor 1 protein (225-296AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-homologous end-joining factor 1 (Protein cernunnos) (XRCC4-like factor), NHEJ1, XLF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NHEJ1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.