



# NHEJ1 Antibody, Biotin conjugated

Product Code	CSB-PA881005LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9Q4
Immunogen	Recombinant Human Non-homologous end-joining factor 1 protein (225-296AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-homologous end-joining factor 1 (Protein cernunnos) (XRCC4-like factor), NHEJ1, XLF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NHEJ1