



NOD2 Antibody, Biotin conjugated

Product Code	CSB-PA881021LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HC29
Immunogen	Recombinant Human Nucleotide-binding oligomerization domain-containing protein 2 protein (460-733AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleotide-binding oligomerization domain-containing protein 2 (Caspase recruitment domain-containing protein 15) (Inflammatory bowel disease protein 1), NOD2, CARD15 IBD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NOD2