



SIRT3 Antibody, Biotin conjugated

Product Code	CSB-PA882102LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTG7
Immunogen	Recombinant Human NAD-dependent protein deacetylase sirtuin-3, mitochondrial protein (20-130AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-dependent protein deacetylase sirtuin-3, mitochondrial (hSIRT3) (EC 3.5.1.-) (Regulatory protein SIR2 homolog 3) (SIR2-like protein 3), SIRT3, SIR2L3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SIRT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.