



MAT2B Antibody, Biotin conjugated

Product Code	CSB-PA882157LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZL9
Immunogen	Recombinant Human Methionine adenosyltransferase 2 subunit beta protein (85-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methionine adenosyltransferase 2 subunit beta (Methionine adenosyltransferase II beta) (MAT II beta) (Putative dTDP-4-keto-6-deoxy-D-glucose 4-reductase), MAT2B, TGR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAT2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.