



STARD9 Antibody, Biotin conjugated

Product Code	CSB-PA885805LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2P6
Immunogen	Recombinant Human StAR-related lipid transfer protein 9 protein (1960-2180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	StAR-related lipid transfer protein 9 (START domain-containing protein 9) (StARD9), STARD9, KIAA1300
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	STARD9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.