



GULP1 Antibody, FITC conjugated

Product Code	CSB-PA887007LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBP9
Immunogen	Recombinant Human PTB domain-containing engulfment adapter protein 1 protein (199-279AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PTB domain-containing engulfment adapter protein 1 (Cell death protein 6 homolog) (PTB domain adapter protein CED-6) (Protein GULP), GULP1, CED6 GULP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GULP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.