




ISCA1 Antibody

Product Code	CSB-PA887126LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUE6
Immunogen	Recombinant Human Iron-sulfur cluster assembly 1 homolog, mitochondrial protein (67-122AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iron-sulfur cluster assembly 1 homolog, mitochondrial (HESB-like domain-containing protein 2) (Iron-sulfur assembly protein IscA) (hIscA), ISCA1, HBLD2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ISCA1

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA887126LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.