



NEUROG2 Antibody, HRP conjugated

Product Code	CSB-PA887960LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2A3
Immunogen	Recombinant Human Neurogenin-2 protein (194-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenin-2 (NGN-2) (Class A basic helix-loop-helix protein 8) (bHLHa8) (Protein atonal homolog 4), NEUROG2, ATOH4 BHLHA8 NGN2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NEUROG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.