



# NEUROG2 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA887960LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9H2A3
<b>Immunogen</b>	Recombinant Human Neurogenin-2 protein (194-262AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Neurogenin-2 (NGN-2) (Class A basic helix-loop-helix protein 8) (bHLHa8) (Protein atonal homolog 4), NEUROG2, ATOH4 BHLHA8 NGN2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	NEUROG2