



NEUROG2 Antibody, Biotin conjugated

Product Code	CSB-PA887960LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2A3
Immunogen	Recombinant Human Neurogenin-2 protein (194-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenin-2 (NGN-2) (Class A basic helix-loop-helix protein 8) (bHLHa8) (Protein atonal homolog 4), NEUROG2, ATOH4 BHLHA8 NGN2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NEUROG2