



WNK1 Antibody, HRP conjugated

Product Code	CSB-PA887978LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4A3
Immunogen	Recombinant Human Serine/threonine-protein kinase WNK1 protein (1-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase WNK1 (EC 2.7.11.1) (Erythrocyte 65 kDa protein) (p65) (Kinase deficient protein) (Protein kinase lysine-deficient 1) (Protein kinase with no lysine 1) (hWNK1), WNK1, HSN2 KDP KIAA0344 PRKWNK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WNK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.