



WNK1 Antibody, Biotin conjugated

Product Code	CSB-PA887978LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4A3
Immunogen	Recombinant Human Serine/threonine-protein kinase WNK1 protein (1-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase WNK1 (EC 2.7.11.1) (Erythrocyte 65 kDa protein) (p65) (Kinase deficient protein) (Protein kinase lysine-deficient 1) (Protein kinase with no lysine 1) (hWNK1), WNK1, HSN2 KDP KIAA0344 PRKWNK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WNK1