



Tropomyosin Antibody, Biotin conjugated

Product Code	CSB-PA889027HD01DRR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9N2R3
Immunogen	Recombinant Charybdis feriata Tropomyosin protein (1-264AA)
Raised In	Rabbit
Species Reactivity	Charybdis feriata
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tropomyosin (Allergen Cha f I) (allergen Cha f 1) (Fragment)
Immunogen Species	Charybdis feriata (Crucifix crab) (Cancer feriatum)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.