



PDS5B Antibody, HRP conjugated

Product Code	CSB-PA889108LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTI5
Immunogen	Recombinant Human Sister chromatid cohesion protein PDS5 homolog B protein (3-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sister chromatid cohesion protein PDS5 homolog B (Androgen-induced proliferation inhibitor) (Androgen-induced prostate proliferative shutoff-associated protein AS3), PDS5B, APRIN AS3 KIAA0979
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PDS5B