



PDS5B Antibody, Biotin conjugated

Product Code	CSB-PA889108LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTI5
Immunogen	Recombinant Human Sister chromatid cohesion protein PDS5 homolog B protein (3-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sister chromatid cohesion protein PDS5 homolog B (Androgen-induced proliferation inhibitor) (Androgen-induced prostate proliferative shutoff-associated protein AS3), PDS5B, APRIN AS3 KIAA0979
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PDS5B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.