



ERVW-1 Antibody, HRP conjugated

Product Code	CSB-PA891578LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQF0
Immunogen	Recombinant Human Syncytin-1 protein (21-443AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syncytin-1 (Endogenous retrovirus group W member 1) (Env-W) (Envelope polyprotein gPr73) (Enverin) (HERV-7q Envelope protein) (HERV-W envelope protein) (HERV-W_7q21.2 provirus ancestral Env polyprotein) (Syncytin) [Cleaved into: Surface protein (SU) (gp50); Transmembrane protein (TM) (gp24)], ERVW-1, ERVWE1
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	ERVW-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.