



# ERVW-1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA891578LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UQF0
<b>Immunogen</b>	Recombinant Human Syncytin-1 protein (21-443AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Syncytin-1 (Endogenous retrovirus group W member 1) (Env-W) (Envelope polyprotein gp73) (Enverin) (HERV-7q Envelope protein) (HERV-W envelope protein) (HERV-W_7q21.2 provirus ancestral Env polyprotein) (Syncytin) [Cleaved into: Surface protein (SU) (gp50); Transmembrane protein (TM) (gp24)], ERVW-1, ERVWE1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Developmental Biology
<b>Target Names</b>	ERVW-1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.