



TMCO1 Antibody, HRP conjugated

Product Code	CSB-PA891784LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM00
Immunogen	Recombinant Human Calcium load-activated calcium channel protein (43-79AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium load-activated calcium channel (CLAC channel) (Transmembrane and coiled-coil domain-containing protein 1) (Transmembrane and coiled-coil domains protein 4) (Xenogeneic cross-immune protein PCIA3), TMCO1, TMCC4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMCO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.