



# TMCO1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA891784LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UM00
<b>Immunogen</b>	Recombinant Human Calcium load-activated calcium channel protein (43-79AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Calcium load-activated calcium channel (CLAC channel) (Transmembrane and coiled-coil domain-containing protein 1) (Transmembrane and coiled-coil domains protein 4) (Xenogeneic cross-immune protein PCIA3), TMCO1, TMCC4
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	TMCO1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.