



ICK Antibody, Biotin conjugated

Product Code	CSB-PA891800LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPZ9
Immunogen	Recombinant Human Serine/threonine-protein kinase ICK protein (150-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase ICK (EC 2.7.11.1) (Intestinal cell kinase) (hICK) (Laryngeal cancer kinase 2) (LCK2) (MAK-related kinase) (MRK), ICK, KIAA0936
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ICK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.