



NOB1 Antibody, FITC conjugated

Product Code	CSB-PA892472LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULX3
Immunogen	Recombinant Human RNA-binding protein NOB1 protein (103-229AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein NOB1 (Phosphorylation regulatory protein HP-10) (Protein ART-4), NOB1, ART4 NOB1P PSMD8BP1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NOB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.