



UBQLN1 Antibody, HRP conjugated

Product Code	CSB-PA892477LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UMX0
Immunogen	Recombinant Human Ubiquilin-1 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquilin-1 (Protein linking IAP with cytoskeleton 1) (PLIC-1) (hPLIC-1), UBQLN1, DA41 PLIC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBQLN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.