



STUB1 Antibody, FITC conjugated

Product Code	CSB-PA892480LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNE7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase CHIP protein (1-48AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase CHIP (EC 2.3.2.27) (Antigen NY-CO-7) (CLL-associated antigen KW-8) (Carboxy terminus of Hsp70-interacting protein) (RING-type E3 ubiquitin transferase CHIP) (STIP1 homology and U box-containing protein 1), STUB1, CHIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STUB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.