



SIGLEC7 Antibody, HRP conjugated

Product Code	CSB-PA896485LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y286
Immunogen	Recombinant Human Sialic acid-binding Ig-like lectin 7 protein (377-467AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sialic acid-binding Ig-like lectin 7 (Siglec-7) (Adhesion inhibitory receptor molecule 1) (AIRM-1) (CDw328) (D-siglec) (QA79 membrane protein) (p75) (CD antigen CD328), SIGLEC7, AIRM1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SIGLEC7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.