



KLF2 Antibody, Biotin conjugated

Product Code	CSB-PA896743LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5W3
Immunogen	Recombinant Human Krueppel-like factor 2 protein (1-355AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Krueppel-like factor 2 (Lung krueppel-like factor), KLF2, LKLF
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	KLF2