



ROBO1 Antibody, Biotin conjugated

Product Code	CSB-PA896760LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6N7
Immunogen	Recombinant Human Roundabout homolog 1 protein (29-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Roundabout homolog 1 (Deleted in U twenty twenty) (H-Robo-1), ROBO1, DUTT1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ROBO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.