



TRAF3IP3 Antibody, HRP conjugated

Product Code	CSB-PA896865LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y228
Immunogen	Recombinant Human TRAF3-interacting JNK-activating modulator protein (74-185AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TRAF3-interacting JNK-activating modulator (TRAF3-interacting protein 3), TRAF3IP3, T3JAM
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRAF3IP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.