



VSIG4 Antibody, HRP conjugated

Product Code	CSB-PA896869LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y279
Immunogen	Recombinant Human V-set and immunoglobulin domain-containing protein 4 protein (20-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	V-set and immunoglobulin domain-containing protein 4 (Protein Z39Ig), VSIG4, CR1g Z39IG
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	VSIG4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.