



SLCO1B1 Antibody, HRP conjugated

Product Code	CSB-PA896932LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6L6
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 1B1 protein (426-537AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 1B1 (Liver-specific organic anion transporter 1) (LST-1) (OATP-C) (Sodium-independent organic anion-transporting polypeptide 2) (OATP-2) (Solute carrier family 21 member 6), SLCO1B1, LST1 OATP1B1 OATP2 OATPC SLC21A6
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SLCO1B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.