




ASB4 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA897287LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y574 |
| Immunogen | Recombinant Human Ankyrin repeat and SOCS box protein 4 protein (166-321AA) |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | ELISA, WB; Recommended dilution: WB:1:500-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ankyrin repeat and SOCS box protein 4 (ASB-4), ASB4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | ASB4 |

Image



Western Blot

Positive WB detected in: A549 whole cell lysate, HepG2 whole cell lysate, Mouse skeletal muscle tissue, Rat heart tissue

All lanes: ASB4 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 49, 39 kDa

Observed band size: 49 kDa