




ASB4 Antibody

Product Code	CSB-PA897287LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y574
Immunogen	Recombinant Human Ankyrin repeat and SOCS box protein 4 protein (166-321AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ankyrin repeat and SOCS box protein 4 (ASB-4), ASB4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ASB4

Image



Western Blot

Positive WB detected in: A549 whole cell lysate, HepG2 whole cell lysate, Mouse skeletal muscle tissue, Rat heart tissue

All lanes: ASB4 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 49, 39 kDa

Observed band size: 49 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.