



DRG1 Antibody, Biotin conjugated

Product Code	CSB-PA897524LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y295
Immunogen	Recombinant Human Developmentally-regulated GTP-binding protein 1 protein (173-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Developmentally-regulated GTP-binding protein 1 (DRG-1) (Neural precursor cell expressed developmentally down-regulated protein 3) (NEDD-3), DRG1, NEDD3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DRG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.