



LYVE1 Antibody, HRP conjugated

Product Code	CSB-PA897574EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5Y7
Immunogen	Recombinant Human Lymphatic vessel endothelial hyaluronic acid receptor 1 protein (20-238AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphatic vessel endothelial hyaluronic acid receptor 1 (LYVE-1) (Cell surface retention sequence-binding protein 1) (CRSBP-1) (Extracellular link domain-containing protein 1) (Hyaluronic acid receptor), LYVE1, CRSBP1 HAR XLKD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LYVE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.