



OAS3 Antibody, Biotin conjugated

Product Code	CSB-PA897591LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6K5
Immunogen	Recombinant Human 2'-5'-oligoadenylate synthase 3 protein (424-565AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2'-5'-oligoadenylate synthase 3 ((2-5')oligo(A) synthase 3) (2-5A synthase 3) (EC 2.7.7.84) (p100 OAS) (p100OAS), OAS3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	OAS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.