





PGAM5 Antibody

Product Code	CSB-PA839338LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HS1
Immunogen	Recombinant Human Serine/threonine-protein phosphatase PGAM5, mitochondrial protein (30-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase PGAM5, mitochondrial (EC 3.1.3.16) (Bcl-XL-binding protein v68) (Phosphoglycerate mutase family member 5), PGAM5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PGAM5

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA839338LA01HU at dilution of 1: 100



Immunofluorescent analysis of HeLa cells using CSB-PA839338LA01HU at dilution of 1: 100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.